

## SPECIALTY WHEELCHAIRS

# REHAB SHOWER COMMUNE

The versatile E&J Folding Rehab Shower Commode chair has earned the reputation of being the premier self-propelled chair for personal hygiene care. This chair features an overall narrow width and a corrosion-resistant stainless steel frame. The chair can be positioned over a standard toilet or used as a portable, self-contained commode.

### Features

- Completely removable full length arms for ease of transfer
- Swingaway, detachable footrests permit close approaches to beds, toilets and wheelchairs
- Recessed seat and frame side contours allow easy access for toileting activities while seated
- 24" rear wheels with handrims and solid tires assist with self propelling
- Removable, hygienic plastic pail and lid are included
- Cushioned vinyl seat is sealed against moisture penetration and vented for optimal foam compression
- Non-corrosive composite footplates
- Maximum weight capacity 250 lb **EVENLY DISTRIBUTED** on models with 14 ½" & 17 ½" back height
- Maximum weight capacity 300 lb **EVENLY DISTRIBUTED** on models with 16" & 19" back height
- Chair comes standard with oval opening seat
- Meets Medicare code E0165



### Warranty

Two (2) year limited warranty



**WARNINGS:** Cancer and Reproductive Harm - [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Color variations may occur due to monitor or printer settings.

GF Health Products, Inc. is an  
ISO 13485:2016  
MDSAP Certified Company  
**MDSAP**  
MEDICAL DEVICE SINGLE AUDIT PROGRAM

© 2010 GF Health Products, Inc. All Rights Reserved.  
Everest & Jennings is a registered trademark of GF Health Products, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.grahamfield.com](http://www.grahamfield.com)

**graham field**

GF1000085RevK21



## Specifications:

		Seat Width		Seat Depth		Seat Height		Back Height from Seat		Rear Wheel diameter		overall Depth		overall Height		overall Width		Weight Without Front Riggings		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
12022010	14.5" Low Back	17	43.18	18	45.72	23	58.42	14.5	36.83	24	60.96	40	101.60	37.5	95.25	25	63.50	43	19.50	250	113.40	1/ea
12022110	17.5" High Back	17	43.18	18	45.72	23	58.42	17.5	44.45	24	60.96	40	101.60	40.5	102.87	25	63.50	45	20.41	250	113.40	1/ea
12023010	16" Low Back	20	50.80	19	48.26	23.5	59.69	16	40.64	24	60.96	42.75	108.59	40	101.60	28.25	71.76	48	21.77	300	136.08	1/ea
12023110	19" High Back	20	50.80	19	48.26	23.5	59.69	19	48.26	24	60.96	42.75	108.59	43	109.22	28.25	71.76	50	22.68	300	136.08	1/ea

## Parts & Options

		Meets HCPCS	UOM
90835070J	Oval cutout seat for 17 1/2" and 14 1/2"	A9900 or A9999	1 ea
90835071J	Oval cutout seat for 16" and 19"	A9900 or A9999	1 ea
90H43170J	Commode pail & lid	A9900 or A9999	1 ea
90761752J	Footrest kit stainless steel (pair)	E0175	1 pr

**WARNINGS:** Cancer and Reproductive Harm - [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Color variations may occur due to monitor or printer settings.

GF Health Products, Inc. is an  
ISO 13485:2016  
MDSAP Certified Company  
**MDSAP**  
MEDICAL DEVICE SINGLE AUDIT PROGRAM

© 2010GF Health Products, Inc. All Rights Reserved.  
Everest & Jennings is a registered trademark of GF Health Products, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.  
[www.grahamfield.com](http://www.grahamfield.com)

**graham field**