

## THREE POSITION RECLINER (574G)


### Features

- The Standard of the Industry and the first recliner specifically designed to improve the comfort and mobility of residents in extended care facilities. The Lumex Three Position Recliner easily adjusts into any one of three positions: sitting, TV or full recline. Recliners offer numerous therapeutic benefits including position changes that can aid in circulation and ease discomfort from respiratory, heart and other ailments.
- A position lock allows only the caregiver to secure the chair in any one of the three easily attained positions. Recliner encourages natural relaxation and improved blood circulation.
- Maximum Weight Capacity: 250 lb EVENLY DISTRIBUTED
- The molded activity table can be fixed in position or left unsecured for easy exit from the recliner.
- Ergonomically designed seat and back with varied density foam construction provides correct support. Articulating headrest quickly adjusts to maximize individual comfort requirements.
- For safer, easier transport, the Lumex Three Position Recliner features two, four inch swivel casters with wheel locks and five inch fixed front wheels.
- Durable welded steel epoxy coated frame is easy to clean and maintain.

### Warranty

- Welded Frame: Lifetime Warranty (Defined as Ten (10) Years)
- Mechanical Components: Three (3) Years
- Upholstered Components: One (1) Year
- Casters: One (1) Year

\*Subject to Warranty Statement Limitations



**WARNINGS:** Cancer and Reproductive Harm - [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Color variations may occur due to monitor or printer settings.



GF Health Products, Inc. is an  
ISO 13485:2016  
MDSAP Certified Company  
**MDSAP**  
(MEDICAL DEVICE SINGLE AUDIT PROGRAM)

© 2009, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.



**graham**  
**field**

GF0900049RevC21