



ROLLATORS

3-WHEEL CRUISER

The GF 3-Wheel Cruiser is designed to provide convenient mobility for those individuals that do not require the extra stability of a 4-wheel rollator. This new light weight aluminum walking aid allows the user to have more freedom in movement to perform daily activities.

Features

- Light weight aluminum frame
- Available in 3 new attractive colors: Pewter Metallic, Metallic Burgundy and Bondi Blue
- Ergonomic hand grips with locking loop brakes
- Adjustable handle height for broad range of users
- 8" black wheels for indoor and outdoor use
- Comes complete with zippered carry pouch
- Instructions included
- Designed to meet the requirements of HCPCS code: E0143
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		Width		Width		overall		minimum		maximum		user		Weight		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg	
609201B	Bondi Blue	14.5	36.83	22	55.88	28	71.12	33	83.82	37	93.98	5' 6" - 6' 5"	1.68 - 1.96	250	113.40	1 ea
609201M	Metallic Burgundy	14.5	36.83	22	55.88	28	71.12	33	83.82	37	93.98	5' 6" - 6' 5"	1.68 - 1.96	250	113.40	1 ea
609201P	Pewter Metallic	14.5	36.83	22	55.88	28	71.12	33	83.82	37	93.98	5' 6" - 6' 5"	1.68 - 1.96	250	113.40	1 ea
	Carry Pouch					14.17	35.99							22	9.98	1 ea



GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company



© 2007, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

Color variations may occur due to monitor or printer settings.

graham
field

GF070002aRevC21