



NON- FOLDING WALKERS

LIGHTWEIGHT RIGID WALKER

Features

- Non-folding design provides extra rigidity for use where space is not an issue
- Angled legs provide maximum patient support in all directions and help resist tipping
- Angled side brace helps support patients as they stand up and sit down
- Fits easily around standard toilet fixtures
- Durable 1" aluminum tubing provides maximum strength while remaining lightweight
- Adjustable to fit a broad range of user heights
- Accommodates all leg attachments
- Designed to meet the requirements of HCPCS code: E0130
- Maximum Weight Capacity 400 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



	Width		Width		overall		Depth		minimum		maximum		user		Weight		UOM		
	inside handgrips		inside rear legs		Width		front to rear		Height		Height		Height Range		Capacity				
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg			
1925A	17	43.18	19.75	50.17	22.5	57.15	20.5	52.07	31	78.74	38	96.52	5' 3"	6' 2"	1.6	1.88	400	181.44	4/cs



GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company



© 2007, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.



+1.770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

Color variations may occur due to monitor or printer settings.



GF070009aRevC21