



CANES

FOLDING CANE

Features

- Cane folds to an 11 inch length for convenient transport and storage
- Derby-style handle is contoured for comfort
- Soft Grip models provide cushioning for added comfort to reduce hand stress and fatigue
- Easy push-button height adjustment from 32.5" - 36.5" with locking collet for added security
- Attractive black frame finish resists marring and scratching
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Limited lifetime warranty



		minimum Height		maximum Height		folded Length		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	
5950A-4	Black	32.5	82.55	36.5	92.71	11	27.94	250	113.40	4/cs