




CANES

ADJUSTABLE OFFSET CANES
ORTHO-EASE[®] GRIP

Features

- Offset handle places the user's weight directly over the base for maximum balance and control
- Attractive finish resists marring and scratching
- Easy push-button height adjustment with locking collet for added security
- Available in a variety of finishes and lengths
- Includes laminated color operating instructions
- Grip angle is designed to adjust to comfortably fit the user's hand and wrist angle and to minimize stress on the joints
- Ortho-Ease[®] grip is ideal for arthritics and for those needing a wider grip for greater support
- Comes complete with convenient wrist strap
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
6327	Bronze	31	78.74	39	99.06	250	113.4	6/cs
6327-1	Bronze	31	78.74	39	99.06	250	113.4	1 ea
6328	Bronze, Extra Long	34.5	87.63	42.5	107.95	250	113.4	6/cs