



CANES

FASHION
OFFSET CANE

Features

- Available in a variety of attractive finishes
- Offset handle places the user's weight directly over the base for maximum balance and control
- Easy push-button height adjustment with locking collet for added security
- Height adjustment 31"- 40"
- Soft, warm and contoured Nitrile grip for maximum comfort
- Sculpted grip reduces the chance of slippage, even when wet
- Comes complete with convenient wrist strap
- Designed to meet the requirements of HCPCS code: E0100
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Limited lifetime warranty



		minimum Height		maximum Height		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
5950BI	Blue Ice	31	78.74	40	101.6	250	113.4	4/cs



Intertek

GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company



© 2007, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.



+1.770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

Color variations may occur due to monitor or printer settings.



GF070005aRevD21