



FOLDING WALKERS

SINGLE RELEASE FOLDING WALKERS
STANDARD RUBBERIZED PVC GRIP

Features

- Single release folding mechanism is designed to aid users with limited hand dexterity
- Durable 1" aluminum tubing provides maximum strength while remaining lightweight
- Raised H-frame design facilitates over-toilet positioning
- Adjustable height fits a broad range of user heights
- Accommodates most Lumex accessories and leg attachments
- Includes laminated color operating instructions
- Handgrip is contoured and textured to provide a better gripping surface
- Designed to meet the requirements of HCPCS code: E0135
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



Handgrip is contoured and textured to provide a better gripping surface

		Width		Width		overall		Depth		Depth		minimum		maximum		user		Weight		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg	
613070A	Adult	17.25	43.82	19.75	50.17	23.25	59.06	19.5	49.53	4.38	11.13	32.5	82.55	39.5	100.33	5' 4" - 6' 2"	1.63 - 1.88	300	136.08	4/cs