



CRUTCHES

UNIVERSAL ALUMINUM CRUTCHES

Features

- Made of lightweight, anodized aluminum
- Reinforced center tube and footpiece provide added strength
- Push button adjustment makes changing height easy
- Height adjustment in 1" increments
- Handgrip position is readily adjusted with wing nuts
- Underarm and handgrip pads are comfortable and durable
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0114
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



Replacement Parts

		UOM
36100	Replacement Tips, Gray	1 pr
361-GRIP	Replacement Handgrips, Gray	1 pr

		minimum Height		maximum Height		minimum Armpad to Grip		maximum Armpad to Grip		user Height Range		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	feet	m	lb	Kg	
3610LF-8	Adult, Latex-Free	45	114.30	53	134.62	15	38.10	19	48.26	5' 2" – 5' 10"	1.57 – 1.78	300	136.08	8 pr/cs
3611LF-8	Tall, Latex-Free	53	134.62	61.5	156.21	16	40.64	20	50.80	5' 10" – 6' 6"	1.78 – 1.98	300	136.08	8 pr/cs
3612LF-8	Youth, Latex-Free	37	93.98	46	116.84	9	22.86	13	33.02	4' 6" – 5' 2"	1.37 – 1.57	300	136.08	8 pr/cs
3613C-8	Combo Pack – 4 Adult, 2 Tall, 2 Youth													8 pr/cs



Intertek

GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company



© 2008, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

Color variations may occur due to monitor or printer settings.

graham
field

GF0800159Rev/E21