



FOLDING WALKERS

colorSelect[®] ADULT WALKER WITH WHEELS

Features

- Sturdy 1" aluminum tubing provides maximum strength while remaining lightweight
- Dual-release folding mechanism enables user to fold walker sides independently
- 5" front fixed wheels
- Complete with glide tips that replace rear leg rubber tips to enhance walker maneuverability on all surfaces
- Designed to meet the requirements of HCPCS code: E0143
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- One (1) year limited warranty



		Width		Width		overall		Depth		Depth		minimum		maximum		Weight Capacity		UOM
		inside handgrips		inside rear legs		Width		front to rear		folded		Height		Height		EVENLY DISTRIBUTED		
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
716270B-2	Blue	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs
716270BK-2	Black	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs
716270G-2	Gold	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs
716270P-2	Plum	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs
716270PK-2	Pink	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs
716270R-2	Red	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs
716270B-2S	Silver	17	43.18	20	50.80	25.5	64.77	20.25	51.44	4	10.16	32	81.28	39	99.06	300	136.08	2/cs