




DAILY LIVING AIDS

TOENAIL CLIPPER

Features

- Precision-ground and nickel-chrome plated
- Includes: file/cleaning tip and thumb pad lever
- Side ground and heat-treated for durable, sharper cutting edges
- Position lugs prevent from slipping



		UOM
1790-2	Toenail Clipper	6/bx



GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company

© 2011, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.



+1.770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

Color variations may occur due to monitor or printer settings.