

HEALTHCARE SEATING

SWING AWAY FOLDING TABLE

Features

- Versatile swing-away table provides convenience to patients
- 360 degree swivel table allows items to be conveniently placed in front of the patient or to their side
- Table swings out of the way for easy patient entry or exit from recliner
- Scratch resistant table top is durable and easy to clean

★ See Chart below for Links, Item Numbers, Specifications, Options and Warranty



Specifications

Specifications																					
Item	Description	Meets HCPCS	Width		Length		Thickness		Height above Seat		Weight Capacity EVENLY DISTRIBUTED		FR566G	FR566DG	FR577RG	FR587W	FR587WD	FR597	FR588W	Made in USA	Warranty in Months
			inch	cm	inch	cm	inch	cm	inch	cm	lb	kg	compatible recliner series								
SAFT10	Swing Away Folding Table	E0184	10.00	25.40	19.00	48.26	0.75	1.905	8.50	21.59	25	1.79	■	■	■	■	■	—	—	■	12
SAFT20													—	—	—	—	—	■	—	■	12
SAFT30													—	—	—	—	—	—	—	■	■
													■ Applies — Does Not Apply								

■ Applies — Does Not Apply