




CLINICAL CARE RECLINERS

POWERED BARIATRIC RECLINER

The Lumex Powered Bariatric Recliner is the perfect solution for larger patients who require a wider recliner with an increased weight capacity. This recliner, with a heavy-gauge steel frame, has a maximum weight capacity of 700 lb when evenly distributed! A four-function hand pendant activates the motor, allowing caregivers to safely recline and reposition patients.

Features

- Quick-release trigger mechanism on push handle allows recliner to swiftly move into Trendelenburg position
- Equipped with rechargeable battery pack for operation if the power goes out
- Four swivel casters facilitate transport and maneuverability in confined areas
- Durable, high impact side tables sit flush with chair arm to conveniently hold prep materials. Side table is smooth and features tapered edges for easy cleaning
- Legrest and backrest raise and lower independently, giving the caregiver complete control over positioning
- Comes complete with head pillow for extra support and patient comfort



Multiple Positions for Enhanced Patient Comfort

Model (xxxx represents fabric color)	overall Width		overall Depth		overall Height		seat Width		seat Depth		seat Height		Arm Height from seat		Width between arms		Width between side tables in upright position		Depth fully reclined (Trendelenburg)		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
FR588Wxxx	38	96.5	59	149.9	54	137.16	26	66.04	24	60.96	27	68.58	8	20.32	24.4	61.98	48	121.92	77	195.6	700	317.51	1 ea



GF Health Products, Inc. is an ISO 13485:2016 MDSAP Certified Company



+1.770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. www.grahamfield.com

Color variations may occur due to monitor or printer settings.

