



## CLINICAL CARE RECLINERS

# ORTHO-BIOTIC<sup>™</sup> II RECLINER

Three reclining positions for patient comfort. Plus a caregiver-activated Trendelenburg (or Shock) Position, is easily actuated from right side by a conveniently located foot pedal.

Legrest-Recliner Lever enables the patient or caregiver to easily operate the legrest and adjust the back angle of the recliner.

Pull Out Patient Foot Plate is easily extended when needed and can be safely stored out of the way when not in use.

Universal Headrest that can easily be adjusted up and down on the back of the recliner.

Four Tente<sup>®</sup> swivel casters facilitate transport and maneuverability in confined areas.

Three total locking casters with one directional locking caster.

Wide selection of fabrics and colors. See the [Swatch Card](#) for details.



Easy access pedal on right side of the recliner allows for activation of the Trendelenburg Position.



Legrest-Recline Lever allows the patient or caregiver to adjust the angle of the back



Footplate is easily extended when needed and can be safely stored out of the way when not in use.

### Options:

- Heat and Massage
- Heat Only feature
- Fully Upholstered Armrests
- Side Table
- IV Pole Holder / IV Pole
- Headrest Cover



Multiple Positions for Enhanced Patient Comfort

Model (xxxx represents fabric colors)	Wood Arm Caps	Fully Upholstered Arms	Heat & Massage	Heat Only	overall Width		overall Depth		overall Height		seat Width		seat Depth		seat Height		Arm Height from Seat		Width between arms		Width between side tables in upright position		Depth fully reclined (Trendelenburg)		Weight Capacity EVENLY DISTRIBUTED		UOM
					inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	
FR597G	●				33	83.8	36.5	92.7	45.5	115.6	23	58.4	20.5	52.1	21	53.3	7.75	19.7	24	61	45	114.3	74	188	350	158.8	1 ea
FR597GH	●		●		33	83.8	36.5	92.7	45.5	115.6	23	58.4	20.5	52.1	21	53.3	7.75	19.7	24	61	45	114.3	74	188	350	158.8	1 ea
FR597GHO	●			●	33	83.8	36.5	92.7	45.5	115.6	23	58.4	20.5	52.1	21	53.3	7.75	19.7	24	61	45	114.3	74	188	350	158.8	1 ea
FR597PxxxxPA		●			33	83.8	36.5	92.7	45.5	115.6	23	58.4	20.5	52.1	21	53.3	7.75	19.7	24	61	45	114.3	74	188	350	158.8	1 ea
FR597PHxxxxPA		●	●		33	83.8	36.5	92.7	45.5	115.6	23	58.4	20.5	52.1	21	53.3	7.75	19.7	24	61	45	114.3	74	188	350	158.8	1 ea
FR597PHOxxxxPA		●		●	33	83.8	36.5	92.7	45.5	115.6	23	58.4	20.5	52.1	21	53.3	7.75	19.7	24	61	45	114.3	74	188	350	158.8	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2014, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.