



HOMECARE/PATIENT ROOM SEATING

EVERYDAY HIP CHAIR

The Lumex Everyday Hip Chair is ideal for people who have arthritis or for patients who have undergone hip or knee surgery. The 27" high seat helps to reduce stress and fatigue on joints while sitting and when entering and exiting the chair.

Features

- High density foam in the upholstered seat and back makes sitting more comfortable for longer periods of time
- New adjustable height footrest provides greater comfort while sitting. Footrest height is easily and securely adjusted by two spring-loaded adjustment knobs located on either side of frame
- Padded full-length arms make sitting and standing easier
- The anti-slip strip on the adjustable footrest provides added safety and is an assist during repositioning
- Sturdy metal frame and padded seat and back are easy to clean and maintain
- Maximum Weight Capacity 350 lb - EVENLY DISTRIBUTED



Footrest can be adjusted to three different heights.

	top of footrest to floor		top of seat to top of footrest	
	inch	cm	inch	cm
High	9.25	23.50	18	45.72
Mid	8.25	20.96	19	48.26
Low	7.25	18.42	20	50.80



Blue Ridge (427)



Jade (857)



Rosewood (863)

Model (xxxx represents fabric color)	overall Width		overall Depth		overall Height		seat Width		seat Depth		seat Height		Arm Height from floor		Arm Height from seat		Width between arms		Shipping Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg			
GF4405xxxx	25	63.5	24.5	62.23	40.75	103.51	19	48.26	16	40.64	27.25	69.22	32	81.28	5	12.7	17.75	45.09	52	23.59	350	158.76	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF0900054Rev.02z