



EXTENDED CARE RECLINER

PREFERRED CARE[®] RECLINER WITH DROP ARM


Lumex Preferred Care Recliner plus a unique drop arm mechanism. Self storing, height adjustable drop arms are simple for the caregiver to activate and enables a safer, easier lateral transfer for the resident and caregiver. Padded armrests provide a soft surface for added comfort. Directional wheel lock on right rear caster reduces "crabbing" of recliner and enables the chair to be moved down hallways with less effort.

Standard Features

- Infinite recliner back patented mechanism provides infinite positioning of the back, ranging from 95 degrees to 145 degrees for increased resident comfort and pressure relief. The reclining back can be locked in any position by the caregiver
- A unique suspension system self-adjusts to each person, cradling the resident from head to heel and virtually eliminating the risk of "bottoming out"
- Extra-wide padded armrests and side panels maximize comfort
- Durable epoxy coated welded steel frame is easy to clean and maintain
- Quick-release upholstery allows for easy cleaning, removal and replacements. Vinyl upholstered back is sulfide stain, mildew and oil resistant
- A variety of accessories are available to customize the chair for specific individual needs

Optional Accessories:

- Anti-Tippers
- IV Pole & Mount
- Headrest Pillow
- Head Support
- Body Bolsters
- Adjustable Footrest
- Foot-Drop Assembly
- Activity Tray Table



Blue Ridge (427)

Doe Skin (851)

Jade (857)

Rosewood (863)

Model (xxxx represents fabric color)	Standard Width Seat	Extra-Wide Seat	Drop-Arm	Tilt-In-Space	overall Width		overall Depth		overall Height		seat Width		seat Depth		seat Height		Arm Height from Seat		Width between arms		Depth fully reclined		Weight Capacity EVENLY DISTRIBUTED		UOM
					inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	
FR565Gxxxx	•				28	71.12	36	91.44	48	121.92	19	48.26	19.5	49.53	20	50.8	8	20.32	19	48.26	66	167.64	250	113.4	1 ea
FR565DGxxxx	•	•			28	71.12	36	91.44	48	121.92	19	48.26	19.5	49.53	22	55.88	8	20.32	19.5	49.53	66	167.64	250	113.4	1 ea
FR565TGxxxx	•		•		28	71.12	36	91.44	48	121.92	19	48.26	19.5	49.53	22	55.88	8	20.32	19.5	49.53	66	167.64	250	113.4	1 ea
FR565WGxxxx		•			32	81.28	37	93.98	48	121.92	24	60.96	19.5	49.53	20	50.8	8.5	21.59	24	60.96	66.5	168.91	400	181.44	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field