



EXTENDED CARE RECLINER

DELUXE PREFERRED CARE® RECLINER WIDE

The Lumex FR565WG Patented recliner mechanism allows for almost infinite positioning of recliner

Unique Position Lock can stabilize recliner position for patient comfort

A unique suspension system self-adjusts to each person, cradling the resident from head to heel and virtually eliminating the risk of "bottoming out"

Extra-wide padded armrests and side panels maximize patient comfort

Durable epoxy coated welded steel frame is easy to clean and maintain

Quick-release upholstery allows for easy cleaning, removal and replacements. Vinyl upholstered back is sulfide stain, mildew and oil resistant.

A variety of accessories are available to customize the chair for specific individual needs

Reduced armrest and seat height enables the recliner to fit under dining/ activity tables

Four TENTE® swivel casters facilitate transport and maneuverability in confined areas. Three total locking casters, with easy to locate black locking tab, lock wheel and swivel with a single motion

Directional wheel lock, with high visibility green locking tab on right rear caster, enables the recliner to track straight while moving

Wide Array of Colors and Fabrics.

Optional Accessories:

- Anti-Tippers
- IV Pole & Mount
- Headrest Pillow
- Head Support
- Body Bolsters
- Adjustable Footrest
- Foot-Drop Assembly
- Activity Tray Table



Model (xxxx represents fabric color)	Standard Width Seat	Extra-Wide Seat	Drop-Arm	Tilt-in-Space	overall Width		overall Depth		overall Height		seat Width		seat Depth		seat Height		Arm Height from Seat		Width between arms		Depth fully reclined		Weight Capacity EVENLY DISTRIBUTED		UOM		
					inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm		lb	Kg
FR565Gxxxx	●				28	71.12	36	91.44	48	121.92	19	48.26	19.5	49.53	20	50.8	8	20.32	19	48.26	66	167.64	250	113.4	1 ea		
FR565DGxxxx	●	●			28	71.12	36	91.44	48	121.92	19	48.26	19.5	49.53	22	55.88	8	20.32	19.5	49.53	66	167.64	250	113.4	1 ea		
FR565TGxxxx	●		●		28	71.12	36	91.44	48	121.92	19	48.26	19.5	49.53	22	55.88	8	20.32	19.5	49.53	66	167.64	250	113.4	1 ea		
FR565WGxxxx	●				32	81.28	37	93.98	48	121.92	24	60.96	19.5	49.53	20	50.8	8.5	21.59	24	60.96	66.5	168.91	400	181.44	1 ea		



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF-0900047/Rev/E22