



## BANDAGES

# MEDICOPASTE® BANDAGE

MEDICOPASTE BANDAGE is intended for use as therapy for peripheral circulatory conditions and certain static conditions of the foot and leg (e.g. foot fractures, sprains and strains)

### Features

- Product should not be used to treat acute pulmonary edema, cellulitis, deep vein thrombosis, arterial insufficiency, infected venous ulcers or if sensitivity to any of the ingredients
- Can be used as "flexible cast" in treatment of sprains, strains and fractures
- Ready to use, easy to apply
- Premium 100% cotton bandage coated with zinc oxide; also available with calamine (ferric oxide)
- Made in U.S.A. in ISO 13485 certified manufacturing facility



		overall Width		overall Length		UOM
		inch	cm	yards	m	
1565 3		3	7.62	10	9.1	12 ea/bx
1565 4		4	10.16	10	9.1	12 ea/bx
1565C-3	with calamine	3	7.62	10	9.1	12 ea/bx
1565C-4	with calamine	4	10.16	10	9.1	12 ea/bx



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2007, GF Health Products, Inc. All Rights Reserved.

**graham**  
**field**

GF070019RevF22