



BANDAGES



UNNA'S BOOT DRESSING

UNNA'S BOOT DRESSING is intended for use as therapy for peripheral circulatory conditions and certain static conditions of the foot and leg (e.g. foot fractures, sprains and strains)

Features

- Product should not be used to treat acute pulmonary edema, cellulitis, deep vein thrombosis, arterial insufficiency, infected venous ulcers or if sensitivity to any of the ingredients
- Can be used as "flexible cast" in treatment of sprains, strains and fractures
- Ready to use, easy to apply
- Premium 100% cotton bandage coated with zinc oxide; also available with calamine (ferric oxide)
- Made in U.S.A. in ISO 13485:2016 certified manufacturing facility



		overall Width		overall Length		UOM
		inch	cm	yards	m	
UNNA3		3	7.62	10	9.1	12 ea/bx
UNNA4		4	10.16	10	9.1	12 ea/bx
UNNA-C3	with calamine	3	7.62	10	9.1	12 ea/bx
UNNA-C4	with calamine	4	10.16	10	9.1	12 ea/bx



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.

graham
field

GF1200056RevE22