



BLADES


CONVENTIONAL STERILE SURGICAL BLADES

Feather conventional sterile surgical blades are individually packaged in an easy-to-peel foil pouch with an inner liner for added protection of blade and user. The double honing process ensures sharp cutting edges with superb durability.

Feather, a blade manufacturer of global renown, utilizes high-precision technology and carefully selected materials. Each Feather scalpel is made in Japan and meets strict ISO standards, resulting in exceptionally high quality and consistent strength and sharpness.

Features

- Heat-treated stainless steel blades provide superior cutting quality, strength, and corrosion resistance
- Gamma-irradiated sterile blades are single use
- Color-coded packaging assists with blade identification
- All blades are individually packaged 100 per dispensing box

Item No.	Stainless Steel Blades	Blade Size	Color Code	UOM
2976#10		10	Red	100/box dispenser
2976#11		11	Pink	100/box dispenser
2976#12		12	Orange	100/box dispenser
2976#15		15	Purple	100/box dispenser
2976#20		20	Yellow	100/box dispenser
2976#21		21	Green	100/box dispenser
2976#22		22	Black	100/box dispenser
2976#23		23	Teal	100/box dispenser
2976#24		24	Orange	100/box dispenser



NOTE: Dispose of scalpel in an approved sharps container in accordance with your facility's sharps disposal procedures



GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

© 2021, GF Health Products, Inc. All Rights Reserved.
GF Health Products, Inc. is the US distributor for FEATHER.
FEATHER is a registered trademark of FEATHER Safety Razor Co., Ltd



+1770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

Color variations may occur due to monitor or printer settings.

grahamfield

GF210007/RevC22