



BLADES

CONVENTIONAL SURGICAL BLADE HANDLES

Features

- Stainless steel handles are designed for smooth loading, functional handling ease, and perfect balance
- Three handle options to accommodate different blades
- Handles are individually packaged and sold



Correct Way to Fit Surgical Blades

- Grasp the blade in the jaws of the forceps or needle holder at a 45-90 degree angle
- Do not grip blade with forceps or needle holder across cutting edge of the blade
- While holding on to the spine of the blade just above the blade slot with forceps or needle holder, insert handle into blade track and slide blade back onto handle until it snaps into position

Item No.	Stainless Steel Handles for Conventional Surgical Blade	Fits Blade Sizes	UOM
2977#3		10, 11, 12, 15	1/ea
2977#4		20, 21, 22, 23, 24	1/ea
2977#7		10, 11, 12, 15	1/ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
GF Health Products, Inc. is the US distributor for FEATHER
FEATHER is a registered trademark of FEATHER Safety Razor Co., Ltd

NOTE: Dispose of scalpel in an approved sharps container in accordance with your facility's sharps disposal procedures

graham field

GF2100004Rev/C22