



BLADES

FEATHER SURGICAL BLADE REMOVER

Features

- Palm-sized disposable blade remover
- Holds between 150 and 300 used blades
- Non-sterile
- Easy removal
- Can be used for any size or type of scalpel handles
- Hygienic for the staff against infection
- Dimensions (mm): 58(H) x 79.3(W) x 25(D)
- Weight (g): 21.7
- Blade Insertion Slot (mm): 17 x 5.2
- Material: Polypropylene



	overall Height		overall Width		overall Depth		Blade Insertion Slot (w x h)		Product Weight		UOM
	inch	mm	inch	mm	inch	mm	inch	mm	oz	g	
2990	2.28	58	3.12	79.3	0.98	25	0.67 x 0.2	17 x 5.2	0.77	21.7	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
GF Health Products, Inc. is the US distributor for FEATHER
FEATHER is a registered trademark of FEATHER Safety Razor Co., Ltd

graham
field

GF2100114RevC22