



BLADES

DISPOSABLE BLOOD LANCET

Features

- Designed for single use only
- Sterilized stainless steel
- Microsharp point provides quick and clean puncture
- Precision engineering assures accuracy to 1/100 mm in thickness
- Superior hardness and sharp point
- Straight side spreads blood easily on glass slide
- Individually wrapped and sealed
- Lancet: 1.5 mm (W) x 3.8 mm (L)
- Made in Japan



	overall Width		overall Legth		Accuracy	UOM
	inch	mm	inch	mm		
3554	0.06	1.50	0.15	3.80	.001 mm	200/bx



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2012, GF Health Products, Inc. All Rights Reserved.
GF Health Products, Inc. is the US distributor for FEATHER.
FEATHER is a registered trademark of FEATHER Safety Razor Co., Ltd

graham
field

GF1200079RevC22