


DIAGNOSTIC INSTRUMENTS

CRILE HEMOSTATIC FORCEPS

Features

- Serrations along the entire length of the jaw
- Ratchets used to lock the forceps and hold blood vessels
- Stainless steel
- Professional use only



		Length		UOM
		inch	cm	
2676	Straight	5.5	14	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
GrafcO is a registered trademark of GF Health Products, Inc.

graham
field

GF2100120 RevC22