

EXAM ROOM SUPPLIES

BINOCULAR LOUPE WITH AUXILIARY LENS

Features

- Binocular loupe comes with auxiliary lens for magnification
- Auxiliary lens swings up and down
- Ideal for work that requires varying degrees of magnification
- Lenses are shatter resistant – optical quality acrylic lenses
- Main parts of binocular loupe are made of ABS (Acrylonitrile Butadiene Styrene)
- Made in USA



	Power with auxiliary lens	W.D.* with auxiliary lens		Power without auxiliary lens	W.D.* without auxiliary lens		UOM
		inch	cm		inch	cm	
1204	4.25X	4	10.2	1.75X	14	35.6	1 ea
1205	4.75X	3	7.6	2.25X	8	20.3	1 ea
1206	5.25X	2	5.1	2.75X	6	15.2	1 ea

* W.D.: Working Distance



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2011, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.