

EXAM ROOM SUPPLIES

HANGING
FRACTURE BEDPAN

Features

- Tapered front for easy placement
- Built-in loop handle makes positioning easy
- Plastic guard prevents spills
- Sturdy, disposable polypropylene
- Portable and lightweight
- Individually wrapped
- Latex free
- Maximum Weight Capacity 331 lb - EVENLY DISTRIBUTED
- Maximum volume: 1,000cc



		maximum Volume		Weight Capacity EVENLY DISTRIBUTED		UOM
		oz	cc	lb	Kg	
2302-B	Gold	33.8	1000	331	150.14	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2008, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.