

EXAM ROOM SUPPLIES

STAINLESS STEEL MEASURING BEAKERS

Features

- Graduated in raised lettering, oz/cc
- Comfort handle
- Stainless steel 18/8 grade 304 with 0.8 mm thickness
- Autoclavable



	Volume		UOM
	oz	cc	
3244	16	500	1 ea
3245	32	1000	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2005, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.

graham
field

GF2100189RevC22