

EXAM ROOM SUPPLIES

# PLUNGERLESS SYRINGES

**Features**

- Ideal for home and hospital use alike
- Heavy-walled resistance flint glass, perfectly retempered and strain-free
- Rubber-compound bulb is specially constructed so liquids cannot enter
- Complete unit can be sterilized repeatedly without fear of breakage
- Catheter point
- Individually boxed

**Caution: This product contains natural rubber latex which may cause allergic reactions**



	Volume		UOM
	oz	ml	
3475	0.5	14.8	1 ea
3476	1	29.6	1 ea
3477	2	59.1	1 ea
3479	4	118.3	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2005, GF Health Products, Inc. All Rights Reserved.  
Grafco is a registered trademark of GF Health Products, Inc.