

EXAM ROOM SUPPLIES

UNLABELED  
DRESSING JARS WITH COVERS

Grafco<sup>®</sup> Unlabeled Flint Glass Dressing Jars are excellent for storing and dispensing a variety of supplies, making them ideal for hospitals and clinics.

**Features**

- Transparent flint glass construction provides easy visual of contents inside
- Each jar includes removable overlapping polished aluminum cover
- Cases include six (6) or twelve (12) flint glass jars complete with polished aluminum lids
- Reusable
- Standard jars not autoclavable
- Pyrex<sup>®</sup> jars are sterilizable and autoclavable
- Not made with natural rubber latex
- One (1) year limited warranty
- Other available options:  
Unlabeled Flint Glass Applicator Jars with Covers (Standard and Pyrex<sup>®</sup>, Sold Individually)



		overall Height		overall Diameter		case Weight		UOM
		inch	cm	inch	cm	lb	kg	
3459	Unlabeled Standard Applicator Jar	6.625	16.8	3	7.6	0.7	0.32	1 ea
3460	Unlabeled Standard Dressing Jar	3	7.6	3	7.6	3.9	1.77	12/cs
3460P	Unlabeled Pyrex Dressing Jar	3	7.6	3	7.6	3.9	1.77	12/cs
3461	Unlabeled Standard Dressing Jar	4	10.2	4	10.2	4.4	2.	6/cs
3462	Unlabeled Standard Dressing Jar	5	12.7	5	12.7	5.1	2.31	6/cs
3463	Unlabeled Standard Dressing Jar	6	15.2	6	15.2	18.4	8.35	6/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2016, GF Health Products, Inc. All Rights Reserved.  
Grafco is a registered trademark of GF Health Products, Inc.