

## EXAM ROOM SUPPLIES

# STAINLESS STEEL WASTE RECEPTACLE

### Features

- Polypropylene liner removable
- Metal handle attached to liner
- Stainless steel #410
- No latex



	overall Width		overall Depth		overall Height		Volume		shipping Weight		UOM
	inch	cm	inch	cm	inch	cm	qrt	l	lb	Kg	
8359	10.6	26.9	10.6	26.9	17	43.2	21	19.9	5.5	2.49	1 ea
8360	10.6	26.9	10.6	26.9	25.98	66	32	30.3	7.3	3.31	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

**+1 770.368.4700**

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2010, GF Health Products, Inc. All Rights Reserved.  
Grafco is a registered trademark of GF Health Products, Inc.

**graham**  
**field**

GF1000055RevC22