

## INFECTION CONTROL

# EYE GOGGLES

### Features

- Splash guard applications
- Indirect ventilation
- Protects against flying fragments, objects, large chips, particles, sand or dirt
- PVC frame
- Adjustable elastic webbing headband contains latex
- Meets requirements of ANSI Z87.1

**Caution:** This product contains natural rubber latex, which may cause allergic reactions.



		UOM
9675	One Size Fits All	24/bx



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

**+1 770.368.4700**

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2013, GF Health Products, Inc. All Rights Reserved.  
Grafco is a registered trademark of GF Health Products, Inc.

**graham**  
**field**

GF1300010RevD22