


LABORATORY

# BOTTLE BRUSH

**Features**

- White nylon bristles on twisted galvanized wire handle

Brushes are not autoclavable



	overall Length		bristle Width		bristle Length		UOM
	inch	cm	inch	cm	inch	cm	
3406	11	27.9	1.75	4.4	4.5	11.4	12 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved.  
GrafcO is a registered trademark of GF Health Products, Inc.