


LABORATORY

# SMALL BORE TEST TUBE BRUSH

**Features**

- White bristles with ½" sponge end
- Galvanized twisted wire handle
- Ideal for cleaning small culture tubes, test tubes, sahli tubes and all small bore glassware

Brushes are not autoclavable



	overall Length		bristle Width		bristle Length		UOM
	inch	cm	inch	cm	inch	cm	
3414	6	15.2	0.5	1.3	3	7.6	12 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved.  
Grafco is a registered trademark of GF Health Products, Inc.