

LABORATORY

# CHENILLE WIRE WICKING STEMS

**Features**

- Pipe Cleaner
- Chenille-wrapped, flexible wire
- Ideal for cleaning laboratory glassware, pipets, funnel stems, stopcocks, or small lumen
- Available in 144-packs or in 50 ft rolls
- Proudly made in the USA

Autoclave or any other sterilization  
is **NOT recommended**



	Length		UOM
	inch	cm	
3560	6	15.2	144/pk
3561	12	30.5	144/pk



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved.  
Grafco is a registered trademark of GF Health Products, Inc.