

LABORATORY

HEMATOCRIT CENTRIFUGE



The Grafco Hematocrit Centrifuge is intended to provide laboratory tests for hematocrit values (HCT) in the fields of biochemistry, microbic immunity genetics, and blood separation, as well as general clinic tests.

The Grafco Hematocrit Centrifuge features a digital timer (up to 99 minutes), fixed speed, brake control, and a reading scale. This unit spins loads at a fixed speed of 12,000 RPM and includes a thermostat to regulate temperature. It can hold up to 24 capillary tubes of size 1.2mm (or 1.5mm) x 75mm.

	overall Depth		overall Width		Height lid closed		Height lid open		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	
410E	10.25	26	9.75	24.8	9.5	24.1	~20	~50.8	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2008, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.