

LAMPS

Q-SERIES UV MAGNIFIER HAND HELD WOODS EXAM LAMPS

Features

- Ideal for visual magnification during fluorescence analysis
- Designed to detect trace amounts of fluorescence through its built-in magnifier lens
- Portable, handheld unit
- Magnification: 3x magnification at 8" focal length
- Impact-resistant polycarbonate housing
- Aluminum reflector maximizes irradiance
- Phenolic handle with a friction-lock swivel joint for comfort and flexibility
- UL Listed
- Made in USA, except bulbs



2200



2203



Tube		Typical Peak UV Intensity at 6" (15.2 cm)	overall Length		overall Width		overall Height		UOM
		μW/cm ²	inch	mm	inch	mm	inch	mm	
2200	(2) 4-Watt Integrally Filtered UV-A Tubes	650	7	17.8	5.75	14.6	2.25	5.7	1 ea
2203	(2) 4-Watt Integrally Filtered UV-A Tubes and (2) 4-Watt White Light Tubes	650	7	17.8	5.75	14.6	2	5.1	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.

graham
field

GF-1300023RevB22