

## LAMPS

# E-SERIES HAND HELD UV-A LAMPS

### Features

- Featuring a silver-anodized aluminum housing that makes the lamp both attractive and durable
- Full range of models offers a selection of UV intensity, size and wattage
- Perfect for applications where a compact lightweight lamp is desired



2212



2211



2213



Tube		Typical Peak UV Intensity at 6" (15.2 cm)	overall Length		overall Width		overall Height		UOM
		$\mu\text{W}/\text{cm}^2$	inch	mm	inch	mm	inch	mm	
2211	(1) 4-Watt integrally filtered tube	800	7	17.8	3.25	8.3	2.5	6.4	1 ea
2212	(1) 6-Watt integrally filtered tube	1,280	10.25	26.	3.25	8.3	2.5	6.4	1 ea
2213	(2) 4-Watt integrally filtered tube	620	7	17.8	3.25	8.3	2.5	6.4	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2010, GF Health Products, Inc. All Rights Reserved.  
Grafco is a registered trademark of GF Health Products, Inc.

**graham**  
**field**

GF100013RevC22