



SPHYGMOMANOMETERS

AUTOMATIC WRIST
BLOOD PRESSURE MONITOR

Blood pressure management is one of the most important things you can do to maintain your health and is essential to preventing heart disease and strokes. More doctors are recommending the use of home blood pressure monitors so that patients can regularly, conveniently check their blood pressure.

Lumiscop's Automatic Wrist Blood Pressure Monitor makes home blood pressure management stress free.

Features

- Automatic Inflation and Deflation
- Digital LCD
- Memory Feature - 2 users 90 sets of readings each
- Average value for the latest 3 readings
- Uses 2 AAA batteries (not included)
- Includes a hard storage case
- Irregular Heartbeat Detection



Specifications:

	overall Width		overall Length		overall Depth		UOM
	inch	mm	inch	mm	inch	mm	
1143	2.6	66	2.75	70	1.3	33	1/ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Lumiscop is a registered trademark of GF Health Products, Inc.