

THERMOMETERS

DIGITAL EAR THERMOMETER

**Features**

- Safely measures temperature in seconds
- Memory recalls last 10 readings
- Measures in Fahrenheit and Celsius
- Water resistant tip- no probe covers necessary
- Color Code Temp Alert Technology
- Temperatures taken are equivalent to those taken orally
- Ideal for any age patient

**Includes**

- Thermometer
- Battery
- Instructions (Spanish and English)



		UOM
2215	Digital Ear Thermometer	1/ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2009, GF Health Products, Inc. All Rights Reserved.  
Lumiscope is a registered trademark of GF Health Products, Inc.