



## THERMOMETERS

# DIGITAL THERMOMETER, °F, FLEXIBLE TIP, 24 BULK PACKAGE

### Features

- Auto Shut Off
- Easy to read LCD Display
- Records and Displays oral, rectal, axillary temperature
- Only reads in Fahrenheit
- Range from 90°F ~ 109.9°F
- Use with a probe cover to help prevent cross infection and contamination



### Specifications:

	Width		Length		Height		Range		Accuracy		Memory	Battery Life	Product Weight		UOM
	inch	cm	inch	cm	inch	cm	F°	C°	F°	C°			oz	g	
1858															
1858CASE	0.7	1.8	4.9	12.4	0.35	0.9	90 - 109.9	32.22 - 43.28	+/- 0.2°, 95.9°-107.6°	+/- 0.11°, 35.5°-42°	Last measured Value	Approximately 200 hours	~.35 including battery	~ 10 including battery	24/bx 240/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2010, GF Health Products, Inc. All Rights Reserved.  
Healthteam is a registered trademark of GF Health Products, Inc.

**graham**  
**field**

GF1000004RevE2z