

HEMOCARE BED ACCESSORIES

# HALF LENGTH CLAMP-ON BED RAIL

The Half Clamp-On Bed Rail provides a full head coverage of 32 inches, designed to offer bed protection and eliminating gaps in between the railing. The non-telescoping tube structures offers a very sturdy rail that can be used as a secondary handle and support.



**Features**

- Improved rail design decreases the spacing between the railing thus reducing the risk of entrapment
- Tool-free installation so it easily mounts to most bed frames using convenient clamp (not for use with "Pan-Style" bed surfaces)
- Intended to define the edge of the sleep surface and to help prevent an individual from inadvertently rolling out of bed
- Designed to be used with the Lumex Patriot™ beds that are commonly used as standard-height or low beds
- Heavy-duty steel construction
- Rails are powder-coated with an environmentally-friendly, long-lasting finish and secure gripping surface
- 3 Year Limited Warranty
- Designed to meet the requirements of HCPCS E0305



**Specifications:**

	overall Length		overall Height		product Weight		UOM
	inch	mm	inch	mm	lb	kg	
GF6650B-1	32	813	18	457	16.8	7.6	1 pr



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2020, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.