



LAMPS

GOOSENECK EXAM LAMPS

**Features**

- Maximum of 60 Watt incandescent, 10 Watt LED equivalent, or 13 Watt CFL equivalent bulb
- 115/120V AC. UL listed 3-wire plug with 8 ft (2.4 m) grounded cord
- Flexible 12" (30.5 cm) gooseneck arm maintains its positioning properly over time and adjusts 360 degrees
- Weighted steel base has smooth rounded lip will not rust, cut into wiring, or mar floors
- Chrome-plated pharmacy shade with baked white enamel bowl for maximum illumination
- Solid cast SOK-IT-GUARD™ Safety Lock, slip proof adjustable lock
- Adjustable Height: 34 ½" – 53 ½" (87.6 – 135.9 cm)
- Overall height (Base, Rod, Gooseneck, Shade): 74" (188 cm)
- Available in a range of configurations
- UL approved
- One (1) year limited warranty



		minimum adjustable Height		maximum adjustable Height		overall Height		UOM
		inch	cm	inch	cm	inch	cm	
1697-1	Chrome-plated Base, SOK-IT-GUARD	34.5	87.6	53.5	135.9	74	188	1 ea
1697-1-220	Chrome-plated Base, SOK-IT-GUARD, 220V	34.5	87.6	53.5	135.9	74	188	1 ea
1697-1C	Chrome-plated Base, Clutch-Collar Lock	34.5	87.6	53.5	135.9	74	188	1 ea
1697-1M	Mobile Base, SOK-IT-GUARD	34.5	87.6	53.5	135.9	74	188	1 ea
1697-1M-220	Mobile Base, SOK-IT-GUARD, 220V	34.5	87.6	53.5	135.9	74	188	1 ea
1697-2	Chrome-plated Base, Hospital Grade Plug, SOK-IT-GUARD	34.5	87.6	53.5	135.9	74	188	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2007, GF Health Products, Inc. All Rights Reserved.  
Grafco is a registered trademark of GF Health Products, Inc.