


DAILY LIVING AIDS

TOENAIL CLIPPER

Features

- Precision-ground and nickel-chrome plated
- Includes: file/cleaning tip and thumb pad lever
- Side ground and heat-treated for durable, sharper cutting edges
- Position lugs prevent from slipping



		UOM
1790-2	Toenail Clipper	6/bx



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2011, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.