

DAILY LIVING AIDS

WOOD TRANSFER BOARDS

Designed to aid the transfer from wheelchair, bed, commode, or chair.

Features

- Solid wood construction
- Lacquer-coated finish for smooth transfer
- Tapered ends for easier transfer
- Available with one or two handles



		overall Width		overall Length		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	lb	Kg	
5241-24A-1	1 Handle	8	20.32	24	60.96	440	199.58	1 ea
5241-30A-1	1 Handle	8	20.32	30	76.2	440	199.58	1 ea
5242-24A-1	2 Handles	8	20.32	24	60.96	440	199.58	1 ea
5242-30A-1	2 Handles	8	20.32	30	76.2	440	199.58	1 ea
5242-32A-1	2 Handles	10	25.4	32	81.28	735	333.39	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2017, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF-1700069RevC22