

DAILY LIVING AIDS

INFLATABLE HAIR SHAMPOOER

Features

- Heavy-duty inflatable plastic ring that supports the neck
- Forms a shallow basin for washing hair, dressing head wounds, or irrigating wounds without moving bedridden patients
- Portable, folds flat for easy storage
- Complete with drainage hose



	overall Width		overall Length		UOM
	inch	cm	inch	cm	
5358	23.5	59.69	22.5	57.15	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2011, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF1100045RevC22