

SUNDRIES

UTILITY APRONS

Features

- Features bound edges, neck band, & side tie tapes
- Latex free Vinyl
- Wash with soap and warm water
- Thickness - 4MIL (0.1MM)



	overall Length		overall Width		overall Thickness		UOM
	inch	cm	inch	cm	mil	mm	
3854	36	91.4	27	68.6	4	0.1	1 ea
3855	42	106.7	36	91.4	4	0.1	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Grafco is a registered trademark of GF Health Products, Inc.