

FOOTSTOOLS

HEAVY-DUTY STEP STOOL

Features

- Robust, steel platform designed to support 1000 lb (453 kg)
- Rubber tipped legs
- One (1) year limited warranty



| | overall Width | | overall Depth | | overall Height | | platform Width | | platform Depth | | Product Weight | | Weight Capacity EVENLY DISTRIBUTED | | UOM |
|--------|---------------|-------|---------------|------|----------------|-------|----------------|-------|----------------|-------|----------------|------|---------------------------------------|--------|------|
| | inch | cm | inch | cm | inch | cm | inch | cm | inch | cm | lb | Kg | lb | Kg | |
| 102BAR | 26 | 66.04 | 20 | 50.8 | 8 | 20.32 | 19 | 48.26 | 14.75 | 37.47 | 5.7 | 2.59 | 1000 | 453.59 | 1 ea |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Intensa is a registered trademark of GF Health Products, Inc.

graham
field