

PHLEBOTOMY CHAIR



PHLEBOTOMY CHAIR

Hausted phlebotomy chair is designed with both the patient and phlebotomist in mind. The sensible and affordable design offers a small footprint and modern style while keeping patient comfort and well-being the highest priority.

- Powder coat frame finish - Foghorn
- Fully upholstered padded arm rests
- One flip arm standard
- Height adjustable arm rests



6200 Series Bariatric

6200-xx upholstered bariatric blood draw/phlebotomy chair with padded arm rest. Dual Flip Arm shown on cover

Specifications:

	6200	
	inch	cm
Overall Width:	34.00	86.36
Overall Depth:	29.00	73.66
Overall Height:	37.50	95.25
Seat Width:	24.50	62.23
Seat Depth:	18.50	46.99
Seat Height:	22.00	64.77
Armrest Height: <i>-low adjustment</i>	27.00	62.23
Armrest Height: <i>-high adjustment</i>	36.00	66.04
Back Width:	25.50	64.77
Back Height:	17.00	43.18
Weight Capacity: <small>EVENLY DISTRIBUTED</small>	lb	kg
	800	362.87
Shipping Weight:	55	24.95
Options: <small>available at additional cost</small>	item number	
	None:	6200-xx
Dual Flip Arms:	6200-DA-xx	

— Does Not Apply

Standard Color Options (-xx)

Available in 24 color choices to complement most decors.

WH Super White	DO Doe	PM Palm	SH Shadow	AR Artichoke	TE Teal
NO Northwoods	DB Dusty Blue	PU Plum	WW Wedgewood	BT British Tan	GR Gray
AL Alabaster	ST Saddle Tan	PL Pearl	TL Tea Leaf	MO Moss	OG Ocean Gray
AQ Aqua Green	BU Blue	VI Violet	NA Navy	WG Warm Gray	BL Black

Swatch Card is available from your representative



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2018, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.

Product Hyperlinks Enabled