

## BARIATRIC PATIENT SEATING


# CLINICAL CARE RECLINER

Designed for patients weighing up to 850 lb (385.5 kg) when evenly distributed, the Gendron Bariatric Clinical Care Recliner features swivel and locking casters for maneuverability. High density foam Back features a 42 degree recline and Elevating Legrest support. The recline mechanism is located at the rear of the Recliner for easy caregiver access.

### Areas of Use

- Dialysis
- Oncology
- IV Therapy
- Patient Room
- Bariatric Care

Item No.	Color	Swatch
715xP1117	Grey	
715xP1199	Dark Blue	
715xP180	Brown	
715xP183	Black	
715xP186	Port Wine	



### Specifications

Specifications			Overall Width		Overall Depth		Overall Height		Seat Width		Seat Depth		Seat Height		Back Height		Shipping or Product Weight		Weight Capacity		Seat Tilt	Backrest Tilt	Footrest Tilt	Made in USA	Warranty
Item	UOM	Meets HCPCS	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	from floor	inch	cm	inch	cm	lb	kg	Evenly Distributed	degrees				
7150P1xxx	1/ea	—	38.25	97.16					29.25	74.30	19.25	48.90						850	385.55		42				36
7155P1xxx	1/ea	—	34.25	87.00					25.25	64.14	19.25	48.90						850	385.55		42				36

■ Standard

■ Optional

■ Available

■ Not Available

■ Accessory

— Does Not Apply

🔗 Active Hyperlink



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

Product Hyperlinks Enabled  
+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.gendroninc.com](http://www.gendroninc.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2020, GF Health Products, Inc. ("Graham-Field") All Rights Reserved.  
Gendron is a registered trademark of GF Health Products, Inc.