



PATIENT TRANSPORT



# GENERAL DUTY STRETCHER

The Gendron General Duty Transport Stretcher is designed for use in the hospital or clinic. Proudly made in the USA, the General Duty G1050 provides secure and comfortable transport throughout the facility. Gendron stretchers are designed with features that provide effortless and reliable patient handling.

## Features

- One-function articulation
  - Raise and lower head section
- Carbon steel with carbon steel sheet metal construction
- Two-position fold-down-type side rails
- Heavy-duty 8" casters with side-mounted lever brakes equipped with Steer/Brake functions
- 2" high-density litter pad with low-shear fabric
- Four IV sockets — one at each corner of stretcher
- IV pole storage under stretcher top for optional telescoping IV pole item 10122



## Dimensions

Stretcher Model	Patient Weight Capacity EVENLY DISTRIBUTED		overall Width		overall Length		Stretcher Top (L x W)		Stretcher Top Height		Side Rail Length		Side Rail Height above Stretcher Top		Caster Size		Floor to Base Clearance		Fowler Angle	Stretcher Shipping Weight	
	lb	kg	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	degrees	lb	kg
G1050	350	158.8	28	71.1	77	30.3	75 x 23	190.5 x 58.4	29.5	74.9	57	144.8	11	27.9	8	20.3	10.5	26.7	≤75	166	75.3



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.gendroninc.com](http://www.gendroninc.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2020, GF Health Products, Inc. All Rights Reserved.  
Gendron is a registered trademark of GF Health Products, Inc.